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Christian G. Giske (Sweden, EUCAST Steering Committee and EARS-Net Coordination Group; chairman), Luis Martinez-Martinez (Spain, EUCAST Steering Committee), Rafael Cantón (Spain, chairman of EUCAST), Stefania Stefani (Italy), Robert Skov (Denmark, EUCAST Steering Committee), Youri Glupczynski (Belgium), Patrice Nordmann (France), Mandy Wootton (UK), Vivi Miriagou (Greece), Gunnar Skov Simonsen (Norway, EARS-Net Coordination Group), Helena Zemlickova (Czech republic, EARS-Net Coordination Group), James Cohen-Stuart (The Netherlands) and Marek Gniadkowski (Poland).

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Christian G. Giske

Chairman of EUCAST- Past Chairman of the subcommittee on detection of resistance mechanisms

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2.5 Referências

1. Canton R, Akova M, Carmeli Y, Giske CG, Glupczynski Y, et al. Rapid evolution and spread of carbapenemases among Enterobacteriaceae in Europe. Clin Microbiol Infect. 2012;18:413-31.
2. Vatopoulos A. High rates of metallo- β -lactamase-producing *Klebsiella pneumoniae* in Greece – a review of the current evidence. Euro Surveill. 2008;13(4). doi:pii: 8023.
3. Albiger B, Glasner C, Struelens MJ, Grundmann H, Monnet DL. Carbapenemase-producing Enterobacteriaceae in Europe: assessment by national experts from 38 countries, May 2015. Euro Surveill. 2015;20(45).
4. European Centres for Disease Prevention and Control (ECDC). Antimicrobial resistance surveillance in Europe 2014. Annual report of the European Antimicrobial Resistance Surveillance Network (EARS-Net)
5. Baraniak A, Izdebski R, Fiett J, Gawryszewska I, Bojarska K, et al. NDM-producing Enterobacteriaceae in Poland, 2012-2014: interregional outbreak of *Klebsiella pneumoniae* ST11 and sporadic cases. J Antimicrob Chemother. 2016;71:85-91.
6. Nordmann P, Naas T, Poirel L. Global spread of Carbapenemase-producing Enterobacteriaceae. Emerg Infect Dis. 2011;17(10):1791-8.
7. Souli M, Galani I, Antoniadou A, Papadomichelakis E, Poulakou G et al. An outbreak of infection due to β -lactamase *Klebsiella pneumoniae* carbapenemase 2-producing *K. pneumoniae* in a Greek University Hospital: molecular characterization, epidemiology, and outcomes. Clin Infect Dis 2010;50:364–73.
8. Bratu S, Landman D, Haag R, Recco R, Eramo A, Alam M, Quale J. Rapid spread of carbapenem-resistant *Klebsiella pneumoniae* in New York City: a new threat to our antibiotic armamentarium. Arch Intern Med. 2005; 165:1430-5.
9. Marchaim D, Navon-Venezia S, Schwaber MJ, Carmeli Y. Isolation of imipenem-resistant *Enterobacter* species: emergence of KPC-2 carbapenemase, molecular characterization, epidemiology, and outcomes. Antimicrob Agents Chemother. 2008;52:1413-8.
10. Queenan AM, Bush K. Carbapenemases: the versatile β -lactamases. Clin Microbiol Rev 2007;20:440–58.
11. Falcone M, Mezzatesta ML, Perilli M, Forcella C, Giordano A et al. Infections with VIM-1 metallo β -lactamase producing *Enterobacter cloacae* and their correlation with clinical outcome. J Clin Microbiol 2009;47: 3514–9.

12. Doumith M, Ellington MJ, Livermore DM, Woodford N. Molecular mechanisms disrupting porin expression in ertapenem-resistant *Klebsiella* and *Enterobacter* spp. clinical isolates from the UK. *J Antimicrob Chemother.* 2009;63:659-67.
13. Yamachika S, Sugihara C, Kamai Y, Yamashita M. Correlation between penicillin-binding protein 2 mutations and carbapenem resistance in *Escherichia coli*. *J Med Microbiol.* 2013;62:429-36.
14. Nordmann P, Gniadkowski M, Giske CG, Poirel L, Woodford N, Miriagou V. Identification and screening of carbapenemase-producing Enterobacteriaceae. *Clin Microbiol Infect.* 2012;18:432-8.
15. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Website with MIC-distributions. (<http://mic.eucast.org/>, accessed on 23 December 2012).
16. Glupczynski Y, Huang TD, Bouchahrouf W, Rezende de Castro R, Bauraing C et al. Rapid emergence and spread of OXA-48-producing carbapenem-resistant Enterobacteriaceae isolates in Belgian hospitals. *Int J Antimicrob Agents.* 2012;39:168-72.
17. Tato M, Coque TM, Ruiz-Garbajosa P, Pintado V, Cobo J et al. Complex clonal and plasmid epidemiology in the first outbreak of Enterobacteriaceae infection involving VIM-1 metallo- β -lactamase in Spain: toward endemicity? *Clin Infect Dis.* 2007;45(9):1171-8.
18. Vading M, Samuelsen O, Haldorsen B, Sundsfjord AS, Giske CG. Comparison of disk diffusion, Etest® and VITEK2 for detection of carbapenemase-producing *Klebsiella pneumoniae* with the EUCAST and CLSI breakpoint systems. *Clin Microbiol Infect.* 2011;17:668-74.
19. Giske CG, Gezelius L, Samuelsen O, Warner M, Sundsfjord A, Woodford N. A sensitive and specific phenotypic assay for detection of metallo- β -lactamases and KPC in *Klebsiella pneumoniae* with the use of meropenem disks supplemented with aminophenylboronic acid, dipicolinic acid and cloxacillin. *Clin Microbiol Infect.* 2011;17:552-6.
20. Doyle D, Peirano G, Lascols C, Lloyd T, Church DL, Pitout JD. Laboratory detection of Enterobacteriaceae that produce carbapenemases. *J Clin Microbiol.* 2012;50:3877-80.
21. Porres-Osante N, de Champs C, Dupont H, Torres C, Mammeri H. Use of avibactam to detect Ambler class A carbapenemases and OXA-48 β -lactamases in Enterobacteriaceae. *Diagn Microbiol Infect Dis.* 2014;79:399-400.
22. Huang TD, Berhin C, Bogaerts P, Glupczynski Y. Evaluation of avibactam-supplemented combination disk tests for the detection of OXA-48 carbapenemase-producing Enterobacteriaceae. *Diagn Microbiol Infect Dis.* 2014;79:252-4.
23. van Dijk K, Voets G, Scharringa J, Voskuil S, Fluit A, Rottier W, Leverstein-Van Hall M, Cohen Stuart J. A disc diffusion assay for detection of class A, B and OXA-48 carbapenemases in Enterobacteriaceae using phenyl boronic acid, dipicolinic acid, and temocillin. *Clin Microbiol Infect* 2013. In press
24. Hartl R, Widhalm S, Kerschner H, Apfalter P. Temocillin and meropenem to discriminate resistance mechanisms leading to decreased carbapenem susceptibility with focus on OXA-48 in Enterobacteriaceae. *Clin Microbiol Infect.* 2013;19:E230-2.
25. Miriagou V, Tzelepi E, Kotsakis SD, Daikos GL, Bou Casals J, Tzouveleakis LS. Combined disc methods for the detection of KPC- and/or VIM-positive *Klebsiella pneumoniae*: improving reliability for the double carbapenemase producers. *Clin Microbiol Infect.* 2013;19:E412-5.
26. Nordmann P, Poirel L, Dortet L. Rapid detection of carbapenemase-producing Enterobacteriaceae. *Emerg Infect Dis.* 2012;18:1503-7.
27. Dortet L, Poirel L, Nordmann P. Rapid Identification of carbapenemase types in Enterobacteriaceae and *Pseudomonas* spp. by using a biochemical test. *Antimicrob Agents Chemother.* 2012;56:6437-40.
28. Vasoo S, Cunningham SA, Kohner PC, Simner PJ, Mandrekar JN, Lolans K, Hayden MK, Patel R. Comparison of a novel, rapid chromogenic biochemical assay, the Carba NP test, with the modified Hodge test for detection of carbapenemase-producing Gram-negative bacilli. *J Clin Microbiol.* 2013;51:3097-101.
29. Tijet N, Boyd D, Patel SN, Mulvey MR, Melano RG. Evaluation of the Carba NP test for rapid detection of carbapenemase-producing Enterobacteriaceae and *Pseudomonas aeruginosa*. *Antimicrob Agents Chemother.* 2013;57:4578-80.
30. Dortet L, Agathine A, Naas T, Cuzon G, Poirel L, Nordmann P. Evaluation of the RAPIDE CR CARBA NP, the Rapid CARB ScreenR and the Carba NP test for biochemical detection of carbapenemase-producing Enterobacteriaceae. *J Antimicrob Chemother.* 2015 ; 70):3014-22.
31. Kabir MH, Meunier D, Hopkins KL, Giske CG, Woodford N. A two-centre evaluation of RAPIDE CR CARBA NP for carbapenemase detection in Enterobacteriaceae, *Pseudomonas aeruginosa* and *Acinetobacter* spp. *J Antimicrob Chemother.* 2016 ;71:1213-6.
32. Huang TD, Berhin C, Bogaerts P, Glupczynski Y. Comparative evaluation of two chromogenic tests for rapid detection of carbapenemase in Enterobacteriaceae and in *Pseudomonas aeruginosa* isolates. *J Clin Microbiol.* 2014;52(8):3060-3.

33. Noel A, Huang TD, Berhin C, Hoebeke M, et al. Comparative Evaluation of Four Phenotypic Tests for Detection of Carbapenemase-Producing Gram-Negative Bacteria. *J Clin Microbiol.* 2017 Feb;55(2):510-518.
34. Pasteran F, Veliz O, Ceriana P, Lucero C, Rapoport M, et al. Evaluation of the Blue-Carba test for rapid detection of carbapenemases in gram-negative bacilli. *J Clin Microbiol.* 2015 ;53(6):1996-8.
35. Compain F, Gallah S, Eckert C, Arlet G, Ramahefasolo A, et al. Assessment of Carbapenem Resistance in Enterobacteriaceae with the Rapid and Easy-to-Use Chromogenic β -Carba Test. *J Clin Microbiol.* 2016;54(12):3065-3068
36. van der Zwaluw K, de Haan A, Pluister GN, Bootsma HJ, de Neeling AJ, et al. The carbapenem inactivation method (CIM), a simple and low-cost alternative for the Carba NP test to assess phenotypic carbapenemase activity in gram-negative rods. *PLoS One.* 2015;10(3):e0123690.
37. Yamada K, Kashiwa M, Arai K, Nagano N, Saito R. Comparison of the Modified-Hodge test, Carba NP test, and carbapenem inactivation method as screening methods for carbapenemase-producing Enterobacteriaceae. *J Microbiol Methods.* 2016;128:48-51.
38. Tijet N, Patel SN, Melano RG. Detection of carbapenemase activity in Enterobacteriaceae: comparison of the carbapenem inactivation method versus the Carba NP test *J Antimicrob Chemother.* 2016 ;71(1):274-6.
39. Hrabak J, Studentova V, Walkova R, Zemlickova H, Jakubu V et al. Detection of NDM-1, VIM-1, KPC, OXA-48, and OXA-162 carbapenemases by matrix-assisted laser desorption ionization-time of flight mass spectrometry. *J Clin Microbiol.* 2012;50:2441-3.
40. Lasserre C, De Saint Martin L, Cuzon G, Bogaerts P, Lamar E, et al. Efficient Detection of Carbapenemase Activity in Enterobacteriaceae by Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry in Less Than 30 Minutes. *J Clin Microbiol.* 2015; 53:2163-71.
41. Papagiannitsis CC, Študentova V, Izdebski R, Oikonomou O, Pfeifer Y, et al. Matrix-assisted laser desorption ionization-time of flight mass spectrometry meropenem hydrolysis assay with NH_4HCO_3 , a reliable tool for direct detection of carbapenemase activity. *J Clin Microbiol.* 2015;53(5):1731.
42. Pasteran F, Denorme L, Ote I, Gomez S, De Belder D, Glupczynski Y, Bogaerts P, Ghiglione B, Power P, Mertens P, Corso A. Rapid identification of OXA-48 and OXA-163 Subfamilies in carbapenem-resistant gram-negative bacilli with a novel immunochromatographic lateral flow assay. *J Clin Microbiol.* 2016;54(11):2832-2836.
43. Dortet L, Jousset A, Sainte-Rose V, Cuzon G, Naas T. Prospective evaluation of the OXA-48 K-SeT assay, an immunochromatographic test for the rapid detection of OXA-48-type carbapenemases. *J Antimicrob Chemother.* 2016;71(7):1834-40.
44. Wareham DW, Shah R, Betts JW, Phee LM, Momin MH. Evaluation of an Immunochromatographic Lateral Flow Assay (OXA-48 K-SeT) for Rapid Detection of OXA-48-Like Carbapenemases in Enterobacteriaceae. *J Clin Microbiol.* 2016;54(2):471-3.
45. Meunier D, Vickers A, Pike R, Hill RL, Woodford N, Hopkins KL. Evaluation of the K-SeT R.E.S.I.S.T. immunochromatographic assay for the rapid detection of KPC and OXA-48-like carbapenemases. *J Antimicrob Chemother.* 2016;71(8):2357-9.
46. Glupczynski Y, Evrard S, Ote I, Mertens P, Huang TD, et al. Evaluation of two new commercial immunochromatographic assays for the rapid detection of OXA-48 and KPC carbapenemases from cultured bacteria. *J Antimicrob Chemother.* 2016;71(5):1217-22.

3. *Enterobacteriaceae* produtoras de β -lactamases de espectro estendido (ESBL)

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3.5 Referências

1. Bush K, Jacoby GA, Medeiros AA. A functional classification scheme for β -lactamases and its correlation with molecular structure. *Antimicrob Agents Chemother.* 1995;39:211-1233.
2. Livermore DM. β -Lactamases in laboratory and clinical resistance. *Clin Microbiol Rev.* 1995;8:557-584.
3. Bradford PA. Extended-spectrum β -lactamases in the 21st century: characterization, epidemiology, and detection of this important resistance threat. *Clin Microbiol Rev.* 2001;14:933-951.
4. Naas T, Poirel L, Nordmann P. 2008. Minor extended-spectrum β -lactamases. *Clin Microbiol Infect.* 2008;14(Suppl1):42-52.
5. Cantón R, Novais A, Valverde A, Machado E, Peixe L, Baquero F, Coque TM. Prevalence and spread of extended-spectrum β -lactamase-producing Enterobacteriaceae in Europe. *Clin Microbiol Infect.* 2008;14(Suppl1):144-153.
6. Livermore DM, Canton R, Gniadkowski M, Nordmann P, Rossolini GM, Arlet G, Ayala J, Coque TM, Kern-Zdanowicz I, Luzzaro F, Poirel L, Woodford N. CTX-M: changing the face of ESBLs in Europe. *J Antimicrob Chemother.* 2007;59:165-174.
7. Carattoli A. Animal reservoirs for extended-spectrum β -lactamase producers. *Clin Microbiol Infect.* 2008;14(Suppl1):117-123.
8. European Centres for Disease Prevention and Control (ECDC). Antimicrobial resistance surveillance in Europe 2014. Annual report of the European Antimicrobial Resistance Surveillance Network (EARS-Net).
9. Gniadkowski M. 2008. Evolution of extended-spectrum β -lactamases by mutation. *Clin Microbiol Infect.* 2008;14(Suppl1):11-32.
10. Livermore DM. Defining an extended-spectrum β -lactamase. *Clin Microbiol Infect.* 2008;14(Suppl1):3-10.
11. European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. EUCAST; 2013. Version 3.0 http://www.eucast.org/clinical_breakpoints/ (last accessed 23 December 2012).
12. Clinical and Laboratory Standards Institute. Performance Standards for Antimicrobial Susceptibility Testing: Twenty-first Informational Supplement M 100-S21. Wayne, PA, USA: CLSI; 2011.
13. Hope R, Potz NA, Warner M, Fagan EJ, Arnold E, Livermore DM. Efficacy of practised screening methods for detection of cephalosporin-resistant Enterobacteriaceae. *J Antimicrob Chemother.* 2007;59(1):110-3.
14. Oliver A, Weigel LM, Rasheed JK, McGowan Jr JE Jr, Raney P, Tenover FC. Mechanisms of decreased susceptibility to cefpodoxime in *Escherichia coli*. *Antimicrob Agents Chemother.* 2002; 46:3829-36.
15. Garrec H, Drieux-Rouzet L, Golmard JL, Jarlier V, Robert J. Comparison of nine phenotypic methods for detection of extended-spectrum β -lactamase production by Enterobacteriaceae. *J Clin Microbiol.* 2011;49:1048-57.
16. Leverstein-van Hall MA, Fluit AC, Paauw A, Box AT, Brisse S, Verhoef J. Evaluation of the EtestR ESBL and the BD Phoenix, VITEK 1, and VITEK 2 automated instruments for detection of extended-spectrum β -lactamases in multiresistant *Escherichia coli* and *Klebsiella* spp. *J Clin Microbiol.* 2002;40:3703-11.
17. Spanu T, Sanguinetti M, Tumbarello M, D'Inzeo T, Fiori B, Posteraro B, Santangelo R, Cauda R, Fadda G. Evaluation of the new VITEK 2 extended-spectrum β -lactamase (ESBL) test for rapid detection of ESBL production in Enterobacteriaceae isolates. *J Clin Microbiol.* 2006;44:3257-62.
18. Thomson KS, Cornish NE, Hong SG, Hemrick K, Herdt C, Moland ES. Comparison of Phoenix and VITEK 2 extended spectrum- β -lactamase detection tests for analysis of *Escherichia coli* and *Klebsiella* isolates with well characterized β -lactamases. *J Clin Microbiol.* 2007;45:2380-4.

19. Drieux L, Brossier F, Sougakoff W, Jarlier V. Phenotypic detection of extended-spectrum β -lactamase production in Enterobacteriaceae: review and bench guide. Clin Microbiol Infect. 2008; 14 (Suppl 1):90-103.
20. Paterson DL, Bonomo RA. Extended-spectrum β -lactamases: a clinical update. Clin Microbiol Rev. 2005;18:657-86.
21. Biedenbach DJ, Toleman M, Walsh TR, Jones RN. Analysis of *Salmonella* spp. with resistance to extended-spectrum cephalosporins and fluoroquinolones isolated in North America and Latin America: report from the SENTRY Antimicrobial Surveillance Program (1997-2004). Diagn Microbiol Infect Dis. 2006;54:13-21.
22. Hirakata Y, Matsuda J, Miyazaki Y, Kamihira S, Kawakami S et al. Regional variation in the prevalence of extended-spectrum β -lactamase-producing clinical isolates in the Asia-Pacific region (SENTRY 1998-2002). Diagn Microbiol Infect Dis. 2005;52:323-9.
23. Jeong SH, Song W, Kim JS, Kim HS, Lee KM. Broth microdilution method to detect extended-spectrum β -lactamases and AmpC β -lactamases in Enterobacteriaceae isolates by use of clavulanic acid and boronic acid as inhibitors. J Clin Microbiol. 2009;47:3409-12.
24. Platteel TN, Cohen Stuart JW, de Neeling AJ, Voets GM, Scharringa J et al. Multi-centre evaluation of a phenotypic extended spectrum β -lactamase detection guideline in the routine setting. Clin Microbiol Infect. 2013;19:70-6.
25. M'Zali FH, Chanawong A, Kerr KG, Birkenhead D, Hawkey PM. Detection of extended-spectrum β -lactamases in members of the family Enterobacteriaceae: comparison of the MAST DD test, the double disc and the EtestR ESBL. J Antimicrob Chemother. 2000;45:881-5.
26. Towne TG, Lewis JS 2nd, Herrera M, Wickes B, Jorgensen JH. Detection of SHV-type extended-spectrum β -lactamase in *Enterobacter* isolates. J Clin Microbiol. 2010;48:298-9.
27. Sturenburg E, Sobottka I, Noor D, Laufs R, Mack D. Evaluation of a new cefepime-clavulanate ESBL EtestR to detect extended-spectrum β -lactamases in an Enterobacteriaceae strain collection. J Antimicrob Chemother. 2004;54:134-8.
28. Nordmann P, Dortet L, Poirel L. Rapid detection of extended-spectrum- β -lactamase-producing Enterobacteriaceae. J Clin Microbiol. 2012;50(9):3016-22.
29. Dortet L, Poirel L, Nordmann P. Rapid detection of ESBL-producing Enterobacteriaceae in blood cultures. Emerg Infect Dis. 2015;21(3):504-7.
30. Renvoise A, Decre D, Amarsy-Guerle R, Huang TD, Jost C, et al. Evaluation of the β -Lacta test, a rapid test detecting resistance to third-generation cephalosporins in clinical strains of Enterobacteriaceae. J Clin Microbiol. 2013;51(12):4012-7.
31. Nukaga M, Mayama K, Crichlow GV, Knox JR. Structure of an extended-spectrum class A β -lactamase from *Proteus vulgaris* K1. J Mol Biol. 2002;317:109-17.
32. Petrella S, Renard M, Ziental-Gelus N, Clermont D, Jarlier V, Sougakoff W. Characterization of the chromosomal class A β -lactamase CKO from *Citrobacter koseri*. FEMS Microbiol Lett. 2006;254:285-92.
33. Sturenburg E, Sobottka I, Noor D, Laufs R, Mack D. Evaluation of a new cefepime-clavulanate ESBL EtestR to detect extended-spectrum β -lactamases in an Enterobacteriaceae strain collection. J Antimicrob Chemother. 2004; 54:134-8.
34. Jacoby GA. AmpC β -lactamases. Clin Microbiol Rev. 2009;22:161-82.
35. Munier GK, Johnson CL, Snyder JW, Moland ES, et al. Positive extended-spectrum- β -lactamase (ESBL) screening results may be due to AmpC β -lactamases more often than to ESBLs. J Clin Microbiol. 2010;48:673-4.
36. Bauernfeind A, Schneider I, Jungwirth R, Sahly H, Ullmann U. A novel type of AmpC β -lactamase, ACC-1, produced by a *Klebsiella pneumoniae* strain causing nosocomial pneumonia. Antimicrob Agents Chemother. 1999;43:1924-31.
37. Polsfuss S, Bloemberg GV, Giger J, Meyer V, Bottger EC, Hombach M. Evaluation of a diagnostic flow chart for detection and confirmation of extended spectrum β -lactamases (ESBL) in Enterobacteriaceae. Clin Microbiol Infect. 2012;18:1194-204.
38. Yang JL, Wang JT, Lauderdale TL, Chang SC. Prevalence of extended-spectrum β -lactamases in *Enterobacter cloacae* in Taiwan and comparison of 3 phenotypic confirmatory methods for detecting extended-spectrum β -lactamase production. J Microbiol Immunol Infect. 2009;42:310-6.
39. Jeong SH, Song W, Kim JS, Kim HS, Lee KM. Broth microdilution method to detect extended-spectrum β -lactamases and AmpC β -lactamases in enterobacteriaceae isolates by use of clavulanic acid and boronic acid as inhibitors. J Clin Microbiol. 2009;47:3409-12.

40. Tsakris A, Poulou A, Themeli-Digalaki K, Voulgari E, Pittaras T, Sofianou D, Pournaras S, Petropoulou D. Use of boronic acid disk tests to detect extended- spectrum β -lactamases in clinical isolates of KPC carbapenemase possessing Enterobacteriaceae. *J Clin Microbiol*. 2009;47:3420-6.
41. March A, Aschbacher R, Dhanji H, Livermore DM, Bottcher A et al. Colonization of residents and staff of a long term- care facility and adjacent acute-care hospital geriatric unit by multiresistant bacteria. *Clin Microbiol Infect*. 2010;16:934-44.
42. Ellington MJ, Ekelund O, Aarestrup FM, et al. The role of whole genome sequencing in antimicrobial susceptibility testing of bacteria: report from the EUCAST Subcommittee. *Clin Microbiol Infect*. 2017;23(1):2-22.

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4.5 Referências

1. Jacoby GA. AmpC beta-lactamases. Clin Microbiol Rev. 2009;22:161-82.
2. Philippon A, Arlet G, Jacoby GA. Plasmid-determined AmpC-type β -lactamases. Antimicrob Agents Chemother. 2002;46:1-11.
3. Beceiro A, Bou G. Class C β -lactamases: an increasing problem worldwide. Rev Med Microbiol. 2004;15:141-152.
4. Empel J, Hrabák J, Kozińska A, Bergerová T, Urbášková P, Kern-Zdanowicz I, Gniadkowski M. DHA1-producing *Klebsiella pneumoniae* in a teaching hospital in the Czech Republic. Microb Drug Resist. 2010;16:291-295.
5. D'Andrea MM, Literacka E, Zioga A, Giani T, Baraniak A, Fiett J, Sadowy E, Tassios PT, Rossolini GM, Gniadkowski M, Miriagou V. Evolution of a multi-drug-resistant *Proteus mirabilis* clone with chromosomal AmpC-type cephalosporinases spreading in Europe. Antimicrob Agents Chemother. 2011;55:2735-2742.
6. Bauernfeind A, Schneider I, Jungwirth R, Sahly H, Ullmann U. A novel type of AmpC β -lactamase, ACC-1, produced by a *Klebsiella pneumoniae* strain causing nosocomial pneumonia. Antimicrob Agents Chemother. 1999;43:1924-31.
7. Edquist P, Ringman M, Liljequist BO, Wisell KT, Giske CG. Phenotypic detection of plasmid-acquired AmpC in *Escherichia coli*-evaluation of screening criteria and performance of two commercial methods for the phenotypic confirmation of AmpC production. Eur J Clin Microbiol Infect Dis. 2013; 32:1205-10.
8. Yagi T, Wachino J, Kurokawa H, Suzuki S, Yamane K et al. Practical methods using boronic acid compounds for identification of class C beta-lactamase-producing *Klebsiella pneumoniae* and *Escherichia coli*. J Clin Microbiol. 2005;43:2551-8.
9. Tenover FC, Emery SL, Spiegel CA, Bradford PA, Eells S et al. Identification of plasmid-mediated AmpC β -lactamases in *Escherichia coli*, *Klebsiella* spp., and *Proteus* spp. can potentially improve reporting of cephalosporin susceptibility testing results. J Clin Microbiol. 2009;47:294-9.
10. Tan TY, Ng LS, He J, Koh TH, Hsu LY. Evaluation of screening methods to detect plasmid-mediated AmpC in *Escherichia coli*, *Klebsiella pneumoniae*, and *Proteus mirabilis*. Antimicrob Agents Chemother. 2009;53:146-9.
11. Martinez-Martinez L. Extended-spectrum β -lactamases and the permeability barrier. Clin Microbiol Infect. 2008;14 Suppl 1:82-9.
12. Halstead FD, Vanstone GL, Balakrishnan I. An evaluation of the Mast D69C AmpC Detection Disc Set for the detection of inducible and derepressed AmpC β -lactamases. J Antimicrob Chemother. 2012;67:2303-4.
13. Ingram PR, Inglis TJ, Vanzetti TR, Henderson BA, Harnett GB, Murray RJ. Comparison of methods for AmpC β -lactamase detection in Enterobacteriaceae. J Med Microbiol. 2011; 60(Pt 6):715-21.
14. Peter-Getzlaff S, Polsfuss S, Poledica M, Hombach M, Giger J et al. Detection of AmpC β -lactamase in *Escherichia coli*: comparison of three phenotypic confirmation assays and genetic analysis. J

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5.1 Referências

1. Giske CG. Contemporary resistance trends and mechanisms for the old antibiotics colistin, temocillin, fosfomicin, mecillinam and nitrofurantoin. *Clin Microbiol Infect.* 2015;21(10):899-905.
2. Olaitan AO, Morand S, Rolain JM. Mechanisms of polymyxin resistance: acquired and intrinsic resistance in bacteria. *Front Microbiol.* 2014; 5:643.
3. Liu YY, Wang Y, Walsh TR, Yi LX, Zhang R, et al. Emergence of plasmid-mediated colistin resistance mechanism MCR-1 in animals and human beings in China: a microbiological and molecular biological study. *Lancet Infect Dis.* 2016 ;16(2):161-8.
4. Skov RL, Monnet DL. Plasmid-mediated colistin resistance (*mcr-1* gene): three months later, the story unfolds. *Euro Surveill.* 2016;21(9).
5. Di Pilato V, Arena F, Tascini C, Cannatelli A, Henrici De Angelis L, et al. *mcr-1.2*, a New *mcr* Variant Carried on a Transferable Plasmid from a Colistin-Resistant KPC Carbapenemase-Producing *Klebsiella pneumoniae* Strain of Sequence Type 512. *Antimicrob Agents Chemother.* 2016; 60:5612-5.
6. Xavier BB, Lammens C, Ruhel R, Kumar-Singh S, Butaye P, Goossens H, Malhotra-Kumar S. Identification of a novel plasmid-mediated colistin-resistance gene, *mcr-2*, in *Escherichia coli*, Belgium, June 2016. *Euro Surveill.* 2016; 21(27).
7. European Centres for Disease Prevention and Control (ECDC). Antimicrobial resistance surveillance in Europe 2014. Annual report of the European Antimicrobial Resistance Surveillance Network (EARS-Net)
8. Hinchliffe P, Yang QE, Portal E, Young T, Li H, et al. Insights into the Mechanistic Basis of Plasmid-Mediated Colistin Resistance from Crystal Structures of the Catalytic Domain of MCR-1. *Sci Rep.* 2017;7: 39392.
9. Bergen PJ, Li J, Rayner CR, Nation RL. Colistin methanesulfonate is an inactive prodrug of colistin against *Pseudomonas aeruginosa*. *Antimicrob Agents Chemother.* 2006;50(6):1953-8.
10. Dafopoulou K, Zarkotou O, Dimitroulia E, Hadjichristodoulou C, Gennimata V, Pournaras S, Tsakris A. Comparative Evaluation of Colistin Susceptibility Testing Methods among Carbapenem-Nonsusceptible *Klebsiella pneumoniae* and *Acinetobacter baumannii* Clinical Isolates. *Antimicrob Agents Chemother.* 2015 Aug;59(8):4625-30.
11. Nordmann P, Jayol A, Poirel L. Rapid Detection of Polymyxin Resistance in Enterobacteriaceae. *Emerg Infect Dis.* 2016;22(6):1038-43.

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6.1 Referências

1. Walsh TR. Emerging carbapenemases: a global perspective. *Int J Antimicrob Agents*. 2010;36 Suppl 3:S8-14.
2. Labarca JA, Salles MJ, Seas C, Guzmán-Blanco M. Carbapenem resistance in *Pseudomonas aeruginosa* and *Acinetobacter baumannii* in the nosocomial setting in Latin America. *Crit Rev Microbiol*. 2016;42:276-92.
3. Higgins PG, Pérez-Llarena FJ, Zander E, Fernández A, Bou G, Seifert H. OXA-235, a novel class D β -lactamase involved in resistance to carbapenems in *Acinetobacter baumannii*. *Antimicrob Agents Chemother*. 2013;57(5):2121-6.
4. Noël A, Huang TD, Berhin C, Hoebeke M, Bouchahrouf W, Yunus S, Bogaerts P, Glupczynski Y. Comparative Evaluation of Four Phenotypic Tests for Detection of Carbapenemase-Producing Gram-Negative Bacteria. *J Clin Microbiol*. 2017;55(2):510-518.
5. Lee K, Lim YS, Yong D, Yum JH, Chong Y. Evaluation of the Hodge test and the imipenem-EDTA double-disk synergy test for differentiating metallo- β -lactamase-producing isolates of *Pseudomonas* spp. and *Acinetobacter* spp. *J Clin Microbiol* 2003; 41, 4623-4629.
6. Hansen F, Hammerum AM, Skov R, Haldorsen B, Sundsfjord A, Samuelsen O. Evaluation of the total MBL confirm kit (ROSCO) for detection of metallo- β -lactamases in *Pseudomonas aeruginosa* and *Acinetobacter baumannii*. *Diagn Microbiol Infect Dis*. 2014;79(4):486-8.
7. Samuelsen O, Buarø L, Giske CG, Simonsen GS, Aasnaes B, Sundsfjord A. Evaluation of phenotypic tests for the detection of metallo- β -lactamase-producing *Pseudomonas aeruginosa* in a low prevalence country. *J Antimicrob Chemother*. 2008;61(4):827-30.
8. Fournier D, Garnier P, Jeannot K, Mille A, Gomez AS, Plésiat P. A convenient method to screen for carbapenemase-producing *Pseudomonas aeruginosa*. *J Clin Microbiol*. 2013;51(11):3846-8.
9. Heinrichs A, Huang TD, Berhin C, Bogaerts P, Glupczynski Y. Evaluation of several phenotypic methods for the detection of carbapenemase-producing *Pseudomonas aeruginosa*. *Eur J Clin Microbiol Infect Dis*. 2015;34(7):1467-74.
10. Kabir MH, Meunier D, Hopkins KL, Giske CG, Woodford N. A two-centre evaluation of RAPIDEC® CARBA NP for carbapenemase detection in Enterobacteriaceae, *Pseudomonas aeruginosa* and *Acinetobacter* spp. *J Antimicrob Chemother*. 2016;71(5):1213-6.

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7.5 Referências

1. Cosgrove SE, Sakoulas G, Perencevich EN, Schwaber MJ, Karchmer AW, Carmeli Y. Comparison of mortality associated with methicillin-resistant and methicillin-susceptible *Staphylococcus aureus* bacteremia: a meta-analysis. Clin Infect Dis. 2003;36:53-9.
2. de Kraker ME, Wolkewitz M, Davey PG, Koller W, Berger J, et al. Clinical impact of antimicrobial resistance in European hospitals: excess mortality and length of hospital stay related to methicillin-resistant *Staphylococcus aureus* bloodstream infections. Antimicrob Agents Chemother. 2011;55:1598-605.
3. de Kraker ME, Davey PG, Grundmann H; BURDEN study group. Mortality and hospital stay associated with resistant *Staphylococcus aureus* and *Escherichia coli* bacteremia: estimating the burden of antibiotic resistance in Europe. PLoS Med. 2011;8(10):e1001104.
4. Chambers HF, Deleo FR. Waves of resistance: *Staphylococcus aureus* in the antibiotic era. Nat Rev Microbiol. 2009;7:629-41.
5. García-Álvarez L, Holden MT, Lindsay H, Webb CR, Brown DF, et al. Methicillin-resistant *Staphylococcus aureus* with a novel *mecA* homologue in human and bovine populations in the UK and Denmark: a descriptive study. Lancet Infect Dis. 2011;11:595-603.
6. Chambers HF. Methicillin resistance in staphylococci: molecular and biochemical basis and clinical implications. Clin Microbiol Rev. 1997;10:781-91.
7. Hososaka Y, Hanaki H, Endo H, et al. Characterization of oxacillin-susceptible *mecA*-positive *Staphylococcus aureus*: a new type of MRSA. J Infect Chemother 2007; 13:79–86.
8. Andrade-Figueiredo M, Leal-Balbino TC. Clonal diversity and epidemiological characteristics of *Staphylococcus aureus*: high prevalence of oxacillin-susceptible *mecA*-positive *Staphylococcus aureus* (OS-MRSA) associated with clinical isolates in Brazil. BMC Microbiol. 2016; 16:115.
9. Proulx MK, Palace SG, Gandra S, Torres B, Weir S, Stiles T, Ellison RT 3rd, Goguen JD. Reversion From Methicillin Susceptibility to Methicillin Resistance in *Staphylococcus aureus* During Treatment of Bacteremia. J Infect Dis. 2016; 213:1041-8.
10. Stegger M, Andersen PS, Kearns A, Pichon B, Holmes MA, Edwards G, Laurent F, Teale C, Skov R, Larsen AR. Rapid detection, differentiation and typing of methicillin-resistant *Staphylococcus aureus* harbouring either *mecA* or the new *mecA* homologue *mecALGA251*. Clin Microbiol Infect. 2012; 4:395-400.
11. Pichon B, Hill R, Laurent F, Larsen AR, Skov RL, Holmes M, Edwards GF, Teale C, Kearns AM. Development of a real-time quadruplex PCR assay for simultaneous detection of *nuc*, Panton-Valentine leucocidin (PVL), *mecA* and homologue *mecALGA251*. J Antimicrob Chemother. 2012; 67:2338-41.

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8.5 Referências

1. Melo-Cristino J, Resina C, Manuel V, Lito L, Ramirez M. First case of infection with vancomycin-resistant *Staphylococcus aureus* in Europe. *Lancet*. 2013;382(9888):205
2. Howden BP, Davies JK, Johnson PDR, Stinear TP, Grayson ML. Reduced vancomycin susceptibility in *Staphylococcus aureus*, including vancomycin-intermediate and heterogeneous vancomycin-intermediate strains: resistance mechanisms, laboratory detection and clinical implications. *Clin Microbiol Rev* 2010;1: 99-139
3. Van Hal SJ, Lodise TP, Paterson DL. The clinical significance of vancomycin minimum inhibitory concentration in *Staphylococcus aureus* infections: a systematic review and meta-analysis. *Clinical Infectious Diseases* 2012; 54: 755-771.
4. Chang HJ, Hsu PC, Yang CC, Siu LK, Kuo AJ, et al. Influence of teicoplanin MICs on treatment outcomes among patients with teicoplanin-treated methicillin resistant *Staphylococcus aureus* bacteraemia: a hospital based retrospective study. *J. Antimicrob Chemother* 2012, 67:736-41.
5. Honda H, Doern CD, Michael-Dunne W Jr, Warren DK. The impact of vancomycin susceptibility on treatment outcomes among patients with methicillin resistant *Staphylococcus aureus* bacteremia. *BMC Infect Dis*. 2011; 5:11:335.
6. Lodise TP, Graves J, Evans A, Graffunder E, Helmecke M, et al. Relationship between vancomycin MIC and failure among patients with methicillin-resistant *Staphylococcus aureus* bacteremia treated with vancomycin. *Antimicrob Agents Chemother* 2008; 52: 3315-20
7. Rojas L, Bunsow E, Munoz P, Cercenado E, Rodriguez-Creixems, Bouza E. Vancomycin MICs do not predict the outcome of methicillin-resistant *Staphylococcus aureus* bloodstream infections in correctly treated patients. *J. Antimicrob Chemother* 2012; 7: 1760-8.
8. Sader HS, Jones RN, Rossi KL, Rybak MJ. Occurrence of vancomycin tolerant and heterogeneous vancomycin resistant strains (hVISA) among *Staphylococcus aureus* causing bloodstream infections in nine USA hospitals. *Antimicrob Chemother*. 2009; 64: 1024-8.
9. Holmes NE, Turnidge JD, Munckhof WJ, Robinson JO, Korman TM, O'Sullivan MV, Anderson TL, Roberts SA, Warren SJ, Gao W, Howden BP, Johnson PD. Vancomycin AUC/MIC ratio and 30-day mortality in patients with *Staphylococcus aureus* bacteremia. *Antimicrob Agents Chemother*. 2013 Apr;57(4):1654-63.
10. Holmes NE, Turnidge JD, Munckhof WJ, Robinson JO, Korman TM, O'Sullivan MV, Anderson TL, Roberts SA, Warren SJ, Gao W, Johnson PD, Howden BP. Vancomycin minimum inhibitory concentration, host comorbidities and mortality in *Staphylococcus aureus* bacteraemia. *Clin Microbiol Infect*. 2013 Dec;19(12):1163-8.
11. Wootton M, Howe RA, Hillman R, Walsh TR, Bennett PM, MacGowan AP. A modified population analysis profile (PAP) method to detect hetero-resistance to vancomycin in *Staphylococcus aureus* in a UK hospital. *J Antimicrob Chemother*. 2001; 47: 399-403
12. Tenover FC, Weigel LM, Appelbaum PC, McDougal LK, Chaitram J, McAllister S, Clark N, Killgore G, O'Hara CM, Jevitt L, Patel JB, Bozdogan B. Vancomycin-resistant

Staphylococcus aureus isolate from a patient in Pennsylvania. Antimicrob Agents Chemother. 2004; 48:275-80

13. Satola SW, Farley MM, Anderson KF, Patel JB. Comparison of Detection Methods for Heteroresistant Vancomycin-Intermediate *Staphylococcus aureus*, with the Population Analysis Profile Method as the Reference Method. J Clin Microbiol. 2011; 49: 177-183.
14. Wootton M, MacGowan AP, Walsh TR, Howe RA. A multicenter study evaluating the current strategies for isolating *Staphylococcus aureus* strains with reduced susceptibility to glycopeptides. J Clin Microbiol. 2007;45:329-32.
15. Voss A, Mouton JW, van Elzakker EP, Hendrix RG, Goessens W, et al. A multi-center blinded study on the efficiency of phenotypic screening methods to detect glycopeptide intermediately susceptible *Staphylococcus aureus* (GISA) and heterogeneous GISA (hGISA). Ann Clin Microbiol Antimicrob. 2007;6:9.

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9.5 Referências

1. Mazuski JE. Vancomycin-resistant enterococcus: risk factors, surveillance, infections, and treatment. *Surg Infect*. 2008;9:567-71.
2. Tenover FC, McDonald LC. Vancomycin-resistant staphylococci and enterococci: epidemiology and control. *Curr Opin Infect Dis*. 2005;18:300-5.
3. Hayden MK, Trenholme GM, Schultz JE, Sahm DF. In vivo development of teicoplanin resistance in a VanB *Enterococcus faecium* isolate. *J Infect Dis*. 1993;167:1224-7.
4. Kawalec M, Gniadkowski M, Kedzierska J, Skotnicki A, Fiett J, Hryniewicz W. Selection of a teicoplanin-resistant *Enterococcus faecium* mutant during an outbreak caused by vancomycin-resistant enterococci with the *vanB* phenotype. *J Clin Microbiol*. 2001;39:4274-82.
5. Holmes NE, Ballard SA, Lam MM, Johnson PD, Grayson ML, Stinear TP, Howden BP. Genomic analysis of teicoplanin resistance emerging during treatment of vanB vancomycin-resistant *Enterococcus faecium* infections in solid organ transplant recipients including donor-derived cases, *J Antimicrob Chemother*. 2013;68(9):2134-9.
6. Depardieu F, Perichon B, Courvalin P. Detection of the *van* alphabet and identification of enterococci and staphylococci at the species level by multiplex PCR. *J Clin Microbiol*. 2004;42:5857-60.
7. Boyd DA, Willey BM, Fawcett D, Gillani N, Mulvey MR. Molecular characterization of *Enterococcus faecalis* N06-0364 with low-level vancomycin resistance harboring a novel D-Ala-D-Ser gene cluster, *vanL*. *Antimicrob Agents Chemother*. 2008;52:2667-72.
8. Xu X, Lin D, Yan G, Ye X, Wu S, Guo Y, Zhu D, Hu F, Zhang Y, Wang F, Jacoby GA, Wang M. *vanM*, a new glycopeptide resistance gene cluster found in *Enterococcus faecium*. *Antimicrob Agents Chemother*. 2010;54:4643-7.
9. Lebreton F, Depardieu F, Bourdon N, Fines-Guyon M, Berger P, Camiade S, Leclercq R, Courvalin P, Cattoir V. DAla-d-Ser VanN-type transferable vancomycin resistance in *Enterococcus faecium*. *Antimicrob Agents Chemother*. 2011;55(10):4606-12.
10. Chen C, Sun J, Guo Y, Lin D, Guo Q, et al. High Prevalence of *vanM* in Vancomycin-Resistant *Enterococcus faecium* Isolates from Shanghai, China. *Antimicrob Agents Chemother*. 2015;59(12):7795-8.
11. Ramotar K, Woods W, Larocque L, Toye B. Comparison of phenotypic methods to identify enterococci intrinsically resistant to vancomycin (VanC VRE). *Diagn Microbiol Infect Dis*. 2000;36:119-24.
12. Sivertsen A, Pedersen T, Larssen KW, Bergh K, Ronning TG, et al. A Silenced *vanA* Gene Cluster on a Transferable Plasmid Caused an Outbreak of Vancomycin-Variable Enterococci. *Antimicrob Agents Chemother*. 2016;60(7):4119-27.
13. Thaker MN, Kalan L, Waglechner N, Eshaghi A, Patel SN, et al. Vancomycin-variable enterococci can give rise to constitutive resistance during antibiotic therapy. *Antimicrob Agents Chemother*.

2015;59(3):1405-10.

14. Grabsch EA, Chua K, Xie S, Byrne J, Ballard SA, Ward PB, Grayson ML. Improved Detection of *vanB2*-Containing *Enterococcus faecium* with Vancomycin Susceptibility by Etest Using Oxgall Supplementation. J. Clin. Microbiol 2008; 46; 1961-4
15. Swenson JM, Clark NC, Sahm DF, Ferraro MJ, Doern G, Hindler J, Jorgensen JH, Pfaller MA, Reller LB, Weinstein MP, et al. Molecular characterization and multilaboratory evaluation of *Enterococcus faecalis* ATCC 51299 for quality control of screening tests for vancomycin and high-level aminoglycoside resistance in enterococci. J Clin Microbiol. 1995;33:3019-21.
16. Klare I, Fleige C, Geringer U, Witte W, Werner G. Performance of three chromogenic VRE screening agars, two EtestR vancomycin protocols, and different microdilution methods in detecting *vanB* genotype *Enterococcus faecium* with varying vancomycin MICs. Diagn Microbiol Infect Dis. 2012;74:171-6.
17. Wijesuriya TM, Perry P, Pryce T, Boehm J, Kay I, Flexman J, Coombs GW, Ingram PR. Low vancomycin MICs and fecal densities reduce the sensitivity of screening methods for vancomycin resistance in Enterococci. J Clin Microbiol. 2014 Aug;52(8):2829-33
18. Endtz HP, Van Den Braak N, Van Belkum A, Goessens WH, Kreft D, Stroebe AB, Verbrugh HA. Comparison of eight methods to detect vancomycin resistance in enterococci. J Clin Microbiol. 1998;36:592-4.
19. Hegstad K, Giske CG, Haldorsen B, Matuschek E, Schonning K, et al. Performance of the EUCAST disk diffusion method, the CLSI agar screen method, and the Vitek 2 automated antimicrobial susceptibility testing system for detection of clinical isolates of Enterococci with low- and medium-level VanB-type vancomycin resistance: a multicenter study. J Clin Microbiol. 2014;52(5):1582-9
20. Fang H, Ohlsson AK, Ullberg M, Ozenci V. Evaluation of species-specific PCR, Bruker MS, VITEK MS and the VITEK2 system for the identification of clinical *Enterococcus* isolates. Eur J Clin Microbiol Infect Dis. 2012; 31: 3073-7
21. Dutka-Malen S, Evers S, Courvalin P. Detection of glycopeptide resistance genotypes and identification to the species level of clinically relevant enterococci by PCR. J Clin Microbiol. 1995;33:1434.
22. Dahl KH, Simonsen GS, Olsvik O, Sundsfjord A. Heterogeneity in the *vanB* gene cluster of genomically diverse clinical strains of vancomycin-resistant enterococci. Antimicrob Agents Chemother. 1999;43:1105-10.
23. Gazin M, Lammens C, Goossens H, Malhotra-Kumar S; MOSAR WP2 Study Team. Evaluation of GeneOhm VanR and Xpert *vanA/vanB* molecular assays for the rapid detection of vancomycin-resistant enterococci. Eur J Clin Microbiol Infect Dis. 2012;31:273-6.

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10.5 Referências

- Kaplan SL, Mason EO Jr. Management of infections due to antibiotic-resistant *Streptococcus pneumoniae*. Clin Microbiol Rev. 1998;11(4):628-44.
- Dagan R. Impact of pneumococcal conjugate vaccine on infections caused by antibiotic-resistant *Streptococcus pneumoniae*. Clin Microbiol Infect. 2009;15 (Suppl 3):16-20.
- Hakenbeck R, Kaminski K, König A, van der Linden M, Paik J, Reichmann P, Zähler D. Penicillin-binding proteins in beta-lactam-resistant *Streptococcus pneumoniae*. Microb Drug Resist 1999; 5: 91-99.
- Grebe T, Hakenbeck R. Penicillin-binding proteins 2b and 2x of *Streptococcus pneumoniae* are primary resistance determinants for different classes of β -lactam antibiotics. Antimicrob Agents Chemother 1996; 40: 829-834.
- Weinstein MP, Klugman KP, Jones RN. Rationale for revised penicillin susceptibility breakpoints versus *Streptococcus pneumoniae*: Coping with antimicrobial susceptibility in an era of resistance. Clin Infect Dis 2009; 48: 1596 – 1600.
- Dixon JMS, Lipinski AE, Graham MEP. Detection and prevalence of pneumococci with increased resistance to penicillin. Can Med Assoc J 1977; 117: 1159-61.
- Swenson JM, Hill BC, Thornsberry C. Screening pneumococci for penicillin resistance. J Clin Microbiol 1986; 24: 749-52.
- Jetté LP and C Sinave. Use of an oxacillin disk screening test for detection of penicillin- and ceftriaxone-resistant pneumococci. J Clin Microbiol 1999; 37: 1178-81